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H_TNFRSF9(4-1BB) Reporter 293 Cell Line

H_TNFRSF9(4-1BB) Reporter 293 细胞系

For research use only!

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Reagent	Specification	Manufacturer/Catalogue No.
Puromycin	25 mg	Genomeditech/GM-040401-1
Blasticidin	10 mg	Genomeditech/GM-040404-1
Pen/Strep	100 mL	Thermo/15140-122
Fetal Bovine Serum	500 mL	ExCell/FSP500
DMEM	500 mL	gibco/C11995500BT
96 Well Clear V-Bottom Tissue Culture	96-well	Corning/3894
96 well round well culture plate	96-well	NEST/701001
96 well White Flat Bottom Polystyrene Not Treated	96-well	Corning/3912
Microplate		
Urelumab (anti-TNFRSF9)	/	Aladdin/Ab170654
Human 4-1BB Ligand/TNFSF9 Trimer Protein	100 µg	kactus/BBL-HM141
Anti-H_4-1BB hIgG2 Antibody(Utomilumab)	/	Genomeditech/GM-26840AB
H_FCGR2B(CD32B) CHO-K1 Cell Line	/	Genomeditech/GM-C16925
GMOne-Step 2.0 Luciferase Reporter Gene Assay Kit	1000T	Genomeditech/GM-040513C

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Equipment	Manufacturer/Catalogue No.
细胞计数仪	ThermoFisher Scientific/Countess 3
酶标仪	Moleculardevices/SpectraMax L

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
Fig 2. Response to Urelumab. The H_TNFRSF9(4-1BB) Reporter 293 Cell line (Cat. GM-C04832) at a concentration of 1.5×10^4 cells/well (96-well format) was stimulated with serial dilutions of Urelumab (anti-TNFRSF9) (Aladdin/Ab170654) in assay buffer (DMEM + 1% FBS + 1% P.S) for 6 hours. The firefly luciferase activity was measured using the Luciferase Reporter Assay Kit (Genomeditech). The maximum induction fold was approximately [3.5].

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
Fig 3. Response to Anti-H_4-1BB hIgG2 Antibody(Utamilumab). H_FCGR2B(CD32B) CHO-K1 Cell Line (Cat. GM-C16925) and CHO-K1 Cell Line were seeded at a density of 1E4 cells/well in a 96-well plate and incubated overnight. The next day, serial dilutions of the Anti-H_4-1BB hIgG2 Antibody(Utamilumab) (Cat. GM-26840AB) were incubated with 1.5E4 cells/well of the H_TNFRSF9(4-1BB) Reporter 293 Cell line (Cat. GM-C04832) in a 96-well plate, and then added to the pre-seeded cells. The mixture was incubated for an additional 6 hours. Firefly

luciferase activity is then measured using the Luciferase Reporter Assay Kit (Genomeditech). The results indicated maximum blocking folds of approximately [10.4]. Data are shown by drug mass concentration.

Genomeditech

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
Fig 4. H_TNFRSF9(4-1BB) Reporter 293 Cell line (Cat. GM-C04832) was determined by flow cytometry using Anti-TNFRSF9(4-1BB) hIgG4 Reference Antibody.

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
Fig 5. The passage stability of the H_TNFRSF9(4-1BB) Reporter 293 Cell line (Cat. GM-C04832) was determined by flow cytometry using Anti-TNFRSF9(4-1BB) hIgG4 Reference Antibody.

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
Fig 6. The passage stability of response to Urelumab. The passages 7, 15, 25, and 32 of H_TNFRSF9(4-1BB) Reporter 293 Cell Line (Cat. GM-C04832) at a density of 1.5E4 cells/well (96-well format) were stimulated with serial dilutions of Urelumab (anti-TNFRSF9) (Aladdin/Ab170654) in assay buffer (DMEM+1% FBS+1% P.S) for 6 hours. The firefly luciferase activity was measured using the Luciferase Reporter Assay Kit (Genomeditech).

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
Fig 7. Response to Human 4-1BB Ligand/TNFSF9 Trimer Protein. The H_TNFRSF9(4-1BB) Reporter 293 Cell line (Cat. GM-C04832) at a density of 1.5E4 cells/well (96-well format) was stimulated with serial dilutions of Human 4-1BB Ligand/TNFSF9 Trimer Protein (kactus/BBL-HM141) in assay buffer (DMEM + 1% FBS + 1% P.S) for 6 hours. The firefly luciferase activity was measured using the Luciferase Reporter Assay Kit

(Genomeditech). The maximum induction fold was approximately [6.0]. Data are shown by drug mass concentration.

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4-1BB	
H_TNFRSF9(4-1BB) Reporter Jurkat Cell Line	H_TNFRSF9(4-1BB) Reporter Jurkat Cell Line (Low Expression)
Cynomolgus_TNFRSF9(4-1BB) CHO-K1 Cell Line	H_TNFRSF9(4-1BB) CHO-K1 Cell Line
Anti-H_4-1BB hIgG2 Antibody(Utomilumab)	
Cynomolgus TNFRSF9(4-1BB) Protein; His Tag	Human TNFRSF9(4-1BB) Protein; His Tag
Human TNFSF9(4-1BBL) Protein; hFc Tag	
CD3	
Cynomolgus_CD3D CD3E Reporter Jurkat(CD3DE KO) Cell Line	Jurkat CD3-BsAb Reporter Cell Line
Cynomolgus_CD3 HEK-293 Cell Line	Cynomolgus_CD3E(Membrane Bound ECD) CHO-K1 Cell Line
H_CD3 CHO-K1 Cell Line	H_CD3 HEK-293 Cell Line
H_CD3(TCR V2) CHO-K1 Cell Line	H_CD3(TCR V2) HEK-293 Cell Line
H_CD3(TCRαβ delta V) HEK-293 Cell Line	H_CD3D CD3E KO Jurkat Cell Line
H_CD3E KO Jurkat Cell Line	H_CD3E(Membrane Bound ECD) CHO-K1 Cell Line
Mouse_CD3 HEK-293 Cell Line	
Anti-CD19×CD3 hIgG1 Antibody[PIT-565(CD58 K34A)]	Anti-CD3 epsilon hIgG1 Antibody [OKT-3 (muromonab)]
Anti-CD3 hIgG1 Antibody(CH2527)	Anti-CD3×CD20 hIgG1 Bispecific Antibody (Epcobio)
Anti-CD3×FCRL5 hIgG1 Bispecific Antibody(cevostamab)	Anti-CD3E×BCMA hIgG4 Reference Antibody (Tecbio)
Anti-CD3E×DLL3 hIgG1 Bispecific Antibody(Tarlabio)	Anti-CD3E×MUC17 hIgG1 Bispecific Antibody(Vepsitbio)
Anti-mouse CD3ε mIgG2a Antibody(145-2C11)	

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