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NFKB Reporter HEK-293 Cell Line

NFKB Reporter HEK-293 细胞系

For research use only!

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Reagent	Specification	Manufacturer/Catalogue No.
Blasticidin	10 mg	Genomeditech/GM-040404-1
Pen/Strep	100 mL	Thermo/15140-122
Fetal Bovine Serum	500 mL	ExCell/FSP500
DMEM	500 mL	Gibco/C11995500BT
96 well round cell culture plate	Plate	PakGent/CL-PT096
96-well White Opaque Plate	Plate	Thermo/236108
96-well U-bottom Plate	Plate	Saining/1014010
GMOne-Step 2.0 Luciferase Reporter Gene Assay Kit	1000T	Genomeditech/GM-040513C
Recombinant Human TNF- α	/	Novoprotein/C008

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Equipment	Manufacturer/Catalogue No.
细胞计数仪	ThermoFisher Scientific/Countess 3
酶标仪	Moleculardevices/SpectraMax L

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

Fig 2. Response to Human TNF- α protein. NFKB Reporter HEK-293 Cell Line (Cat. GM-C29472) at a concentration of 1.5E4 cells/well (96-well format) was stimulated with serial dilutions of Recombinant Human TNF- α (Novoprotein/C008) in assay buffer (DMEM + 1% FBS + 1% P.S) for 7 hours. The firefly luciferase activity was measured using the Luciferase Reporter Assay Kit (Genomeditech). The maximum induction fold was approximately [22.8].

(Right panel: Data plotted after conversion of drug concentrations from mass concentration to molar concentration.)

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
Fig 3. The passage stability of response to Human TNF- α protein. The passages 6, 21, and 36 of NFKB Reporter HEK-293 Cell Line (Cat. GM-C29472) at a density of 1.5E4 cells/well (96-well format) were stimulated with serial dilutions of Recombinant Human TNF- α (Novoprotein/C008) in assay buffer (DMEM + 1% FBS + 1% P.S) for 7 hours. The firefly luciferase activity was measured using the Luciferase Reporter Assay Kit (Genomeditech).

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NFKB	
NFKB Reporter A549 Cell Line	NFKB Reporter HEK-293 Cell Line
NFKB Reporter Jurkat Cell Line	NFKB Reporter TF-1 Cell Line
NFKB Reporter THP1 Cell Line	

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