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### H\_TREM1 Blockade Reporter Jurkat Cell Line

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Reagent	Specification	Manufacturer/Catalogue No.
Puromycin	25 mg	Genomeditech/GM-040401-1
Blasticidin	10 mg	Genomeditech/GM-040404-1
Pen/Strep	100 mL	Thermo/15140-122
Fetal Bovine Serum	500 mL	ExCell/FSP500
RPMI 1640	500 mL	Gibco/C11875500BT
96 Well Clear V-Bottom Tissue Culture	96-well	Corning/3894
Clear Flat-Bottom Immuno Nonsterile 96-Well Plates	96-well	Thermo/442404
96 well White Flat Bottom Polystyrene Not Treated	96-well	Corning/3912
Microplate		
Anti-H_TREM1 hIgG1 Antibody	/	Genomeditech/GM-26835AB
Human PGLYRP1 Protein; His Tag	/	Genomeditech/GM-87969RP
Peptidoglycan	/	ABMole/M21611
H_PGLYRP1 HEK-293 Cell Line	/	Genomeditech/GM-C41231
Anti-TREM1 hIgG1 Antibody(mAb0170)	/	Genomeditech/GM-26834AB
GMOne-Step 2.0 Luciferase Reporter Gene Assay Kit	1000T	Genomeditech/GM-040513C

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Equipment	Manufacturer/Catalogue No.
细胞计数仪	ThermoFisher Scientific/Countess 3
酶标仪	Moleculardevices/SpectraMax L

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
Fig 2. Response to Human PGLYRP1 Protein; His Tag. The H\_TREM1 Blockade Reporter Jurkat Cell Line (Cat. GM-C15720) at a concentration of 1E5 cells/well (96-well format) was stimulated with serial dilutions of Human PGLYRP1 Protein; His Tag (Cat. GM-87969RP) and 1 µg/mL Peptidoglycan (ABMole/M21611) in assay buffer (RPMI 1640 + 1% FBS + 1% P.S) for 6 hours. The firefly luciferase activity was measured using the Luciferase Reporter Assay Kit (Genomeditech). The maximum induction fold was approximately [43.7]. Data are shown by drug mass concentration.

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
Fig 3. The passage stability of response to Human PGLYRP1 Protein; His Tag. The passage 4, 15, and 32 of H\_TREM1 Blockade Reporter Jurkat Cell Line (Cat. GM-C15720) at a concentration of 1E5 cells/well (96-well format) were stimulated with serial dilutions of Human PGLYRP1 Protein; His Tag (Cat. GM-87969RP) and 1 µg/mL Peptidoglycan (ABMole/M21611) in assay buffer (RPMI 1640 + 1% FBS + 1% P.S) for 6 hours. The firefly luciferase activity was measured using the Luciferase Reporter Assay Kit (Genomeditech). Data are shown by drug mass concentration.

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
Fig 4. The passage stability of the H\_TREM1 Blockade Reporter Jurkat Cell Line (Cat. GM-C15720) was determined by flow cytometry using Anti-H\_TREM1 hIgG1 Antibody (Cat. GM-26835AB).

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
Fig 5. Response to Anti-H\_TREM1 hIgG1 Antibody. H\_TREM1 Blockade Reporter Jurkat Cell Line (Cat. GM-C15720) at a concentration of 1E5 cells/well in a 96-well format. The wells were coated overnight with serial dilutions of Anti-H\_TREM1 hIgG1 Antibody (Cat. GM-26835AB) in assay buffer (RPMI 1640 + 1% FBS + 1% P.S). After coating, the cells were added and incubated for 24 hours. The firefly luciferase activity was measured using the Luciferase Reporter Assay Kit (Genomeditech). The maximum induction fold was approximately [178.8]. Data are shown by drug mass concentration.

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
Fig 6. Response to Anti-H\_TREM1 hIgG1 Antibody. Anti-H\_TREM1 hIgG1 Antibody (Cat. GM-26835AB) was seeded at a density of 0.18 µg/well in a 96-well plate and incubated overnight. The next day, serial dilutions of the Anti-H\_TREM1 hIgG1 Antibody (Cat. GM-26835AB) were incubated with 1E5 cells/well of the H\_TREM1 Blockade Reporter Jurkat Cell Line (Cat. GM-C15720) in a 96-well plate, and then added to the pre-seeded plate. The mixture was incubated for an additional 24 hours. Firefly luciferase activity is then measured using the Luciferase Reporter Assay Kit (Genomeditech). The results indicated maximum blocking folds of approximately [9.8]. Data are shown by drug mass concentration.

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
Fig 7. Response to the H\_PGLYRP1 HEK-293 Cell Line. The H\_PGLYRP1 HEK-293 Cell Line (Cat. GM-C41231) at a concentration of 5E4 cells/well was co-cultured with H\_TREM1 Blockade Reporter Jurkat Cell Line (Cat. GM-C15720) at a concentration of 1E5 cells/well, in assay buffer (RPMI 1640 + 1% FBS + 1% P.S) for 6 hours (96-well format). Firefly luciferase activity is then measured using the Luciferase Reporter Assay Kit (Genomeditech). The maximum induction fold was approximately [79.1].

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
Fig 8. Inhibition of H\_PGLYRP1 HEK-293 Cell Line cell-induced reporter activity by Anti-TREM1 hIgG1 Antibody(mAb0170). H\_PGLYRP1 HEK-293 Cell Line (Cat. GM-C41231) was seeded at a density of 1.5E4 cells per well in a 96-well plate and incubated overnight. The next day, serial dilutions of the Anti-TREM1 hIgG1 Antibody(mAb0170) (Cat. GM-26834AB) were incubated with 1E5 cells/well of the H\_TREM1 Blockade Reporter Jurkat Cell Line (Cat. GM-C15720) in a 96-well plate for 1 hour, and then added to the pre-seeded cells. The mixture was incubated for an additional 6 hours. Firefly luciferase activity is then measured using the Luciferase Reporter Assay Kit (Genomeditech)(left Y-axis, relative luminescence units), with inhibition percentages shown on the right Y-axis.

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TREM1	
<a href="#">Cynomolgus_TREM1 CHO-K1 Cell Line</a>	<a href="#">Cynomolgus_TREM1 HEK-293 Cell Line</a>
<a href="#">H_TREM1 CHO-K1 Cell Line</a>	<a href="#">H_TREM1 HEK-293 Cell Line</a>
<a href="#">Mouse_TREM1 CHO-K1 Cell Line</a>	
<a href="#">Anti-TREM1 hIgG1 Antibody</a>	
<a href="#">Human PGLYRP1 Protein; His Tag</a>	
TREM2	
<a href="#">H_TREM2 Reporter Jurkat Cell Line</a>	<a href="#">Cynomolgus_TREM2 CHO-K1 Cell Line</a>
<a href="#">Cynomolgus_TREM2 HEK-293 Cell Line</a>	<a href="#">H_TREM2 CHO-K1 Cell Line</a>
<a href="#">H_TREM2 HEK-293 Cell Line</a>	<a href="#">Mouse_TREM2 HEK-293 Cell Line</a>
<a href="#">Anti-H_TREM2 hIgG4 Antibody</a>	<a href="#">Anti-H_TREM2 Rat_IgG2b Antibody</a>
<a href="#">Anti-TREM2 hIgG1 Antibody</a>	

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