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H_IGF-1R Reporter Cell Line

H_IGF-1R Reporter 细胞系

For research use only!

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Reagent	Specification	Manufacturer/Catalogue No.
DMEM	500 mL	gibco/C11995500BT
Fetal Bovine Serum	500 mL	ExCell/FSP500
Pen/Strep	100 mL	Thermo/15140-122
Blasticidin	10 mg	Genomeditech/GM-040404-1
G418	1 g	Genomeditech/ GM-040402-1
Puromycin	25 mg	Genomeditech/GM-040401-1
96 Well Clear V-Bottom Tissue Culture	96-well	Corning/3894
96 well round well culture plate	96-well	NEST/701001
96 well White Flat Bottom Polystyrene Not Treated Microplate	96-well	Corning/3912
Human IGF-1 Protein; His Tag	/	Genomeditech/GM-87635RP
Anti-IGF-1R hIgG1 Antibody(teprotumumab)	/	Genomeditech/GM-87441AB
GMOne-Step 2.0 Luciferase Reporter Gene Assay Kit	1000T	Genomeditech/GM-040513C

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Equipment	Manufacturer/Catalogue No.
细胞计数仪	ThermoFisher Scientific/Countess 3
酶标仪	Moleculardevices/SpectraMax L

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
Fig 2. Response to Recombinant Human IGF1 Protein. The H_IGF-1R Reporter Cell Line (Cat. GM-C28412) at a concentration of 1.5E4 cells/well (96-well format) was stimulated with serial dilutions of Human IGF-1 Protein; His Tag (Cat. GM-87635RP) in assay buffer (DMEM+1% FBS+1% P.S) for 16 hours. The firefly luciferase

activity was measured using the Luciferase Reporter Assay Kit (Genomeditech). The maximum induction fold was approximately [18.7]. Data are shown by drug mass concentration.

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
Fig 3. Response to Anti-H_IGF-1R hIgG1 Antibody(teprotumumab). Serial dilutions of the Anti-H_IGF-1R hIgG1 Antibody(teprotumumab) were incubated with 1.5E4 cells/well of the H_IGF-1R Reporter Cell Line (Cat.GM-C28412) in a 96-well plate for 1 hour in assay buffer (DMEM+1% FBS+1% P.S). Subsequently, the Recombinant Human IGF1 Protein (Company A) at a concentration of 300 ng/mL was added, and the coculture proceeded for an additional 16 hours. Firefly luciferase activity was then measured using the Luciferase Reporter Assay Kit (Genomeditech). The results indicated maximum blocking folds of approximately [11.5]. Data are shown by drug mass concentration.

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
Fig 4. The passage stability of response to Human IGF1 Protein. The passage 8, 13, 20, 24, and 30 of H_IGF-1R Reporter Cell Line (Cat. GM-C28412) at a concentration of 1.5E4 cells/well (96-well format) were stimulated with serial dilutions of Recombinant Human IGF1 Protein (Company A) in assay buffer (DMEM + 1% FBS + 1% P.S) for 16 hours. The firefly luciferase activity was measured using the Luciferase Reporter Assay Kit (Genomeditech). Data are shown by drug mass concentration.

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
Fig 5. The passage stability of the H_IGF-1R Reporter Cell Line (Cat. GM-C28412) was determined by flow cytometry using Anti-IGF-1R hIgG1 Antibody(teprotumumab) (Cat. GM-87441AB).

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IGF-1R	
H_IGF-1R CHO-K1 Cell Line	H_IGF-1R HEK-293 Cell Line
Anti-IGF-1R hIgG1 Antibody(teprotumumab)	

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