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### H\_TACI Reporter Cell Line

### H\_TACI Reporter 细胞系

For research use only!

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Reagent	Specification	Manufacturer/Catalogue No.
Puromycin	25 mg	Genomeditech/GM-040401-1
Blasticidin	10 mg	Genomeditech/GM-040404-1
Pen/Strep	100 mL	Thermo/15140-122
Fetal Bovine Serum	500 mL	ExCell/FSP500
RPMI 1640	500 mL	Gibco/C11875500BT
96 well round cell culture plate	Plate	PakGent/CL-PT096
96-well White Opaque Plate	Plate	Thermo/236108
96-well U-bottom Plate	Plate	Saining/1014010
H_TACI Reporter Cell Line	1 管(5E6 Cells/mL)	Genomeditech/GM-C35041
Human BAFF Protein; His Tag	/	Genomeditech/GM-87735RP
Anti-BAFF hIgG1 Antibody(belimumab)	/	Genomeditech/GM-87690AB
APC anti-human CD267 (TACI) Antibody	25 tests	BioLegend/311911
Human APRIL / TNFSF13 Protein, His,Flag Tag	100 µg	ACRO/APL-H52D1
Anti-APRIL hIgG2 Reference Antibody (Sibebio)	/	Genomeditech/GM-88014MAB
GMOne-Step 2.0 Luciferase Reporter Gene Assay Kit	1000T	Genomeditech/GM-040513C

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Equipment	Manufacturer/Catalogue No.
细胞计数仪	ThermoFisher Scientific/Countess 3
酶标仪	Moleculardevices/SpectraMax L

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H_TACI Reporter Cell Line	0 $\mu\text{g}/\text{mL}$	5 $\mu\text{g}/\text{mL}$	76.29 $\mu\text{g}/\text{mL}$
		7257	128191

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
Fig 2. Response to Human BAFF Protein. The H\_TACI Reporter Cell Line (Cat. GM-C35041) at a concentration of  $1 \times 10^5$  cells/well (96-well format) was stimulated with serial dilutions of Human BAFF Protein (Cat. GM-87735RP) in assay buffer (RPMI 1640 + 1% FBS + 1% P.S) for 6 hours. The firefly luciferase activity was measured using the Luciferase Reporter Assay Kit (Genomeditech). The maximum induction fold was approximately [18.2].

## 2. Human BAFF Protein 激活; Anti-BAFF 抑制实验

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
Fig 3. Response to Anti-BAFF hlgG1 Antibody(belumumab). Serial dilutions of Anti-BAFF hlgG1 Antibody(belumumab) (Cat. GM-87690AB) was incubated with 4.7 ng/well of Human BAFF Protein (Cat. GM-87735RP) for 1 hour in assay buffer (RPMI 1640 + 1% FBS + 1% P.S). After pre-incubation, add the mixture to the H\_TAC1 Reporter Cell Line (Cat. GM-C35041) at a density of 1E5 cells/well in a 96-well format, and incubate for 6 hours. Firefly luciferase activity is then measured using the Luciferase Reporter Assay Kit (Genomeditech).

The results indicated maximum blocking folds of approximately [6.5].

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
Fig 4. Inhibition of Membrane Bound H\_APRIL-induced reporter activity by Sibeprenlimab. Serial dilutions of the Anti-APRIL hIgG2 Reference Antibody (Sibebio)(Cat. GM-88014MAB) were incubated with 1.5E4 cells/well of the Membrane Bound H\_APRIL(Trimer) HEK-293 Cell Line (Cat. GM-C40503) in a 96-well plate for 1 hour in assay buffer (RPMI 1640+1% FBS+1% P.S). Subsequently, the H\_TACI Reporter Cell Line (Cat. GM-C35041) at a concentration of 1E5 cells/well was added, and the co-culture proceeded for an additional 6 hours. Firefly luciferase activity was then measured using the Luciferase Reporter Assay Kit (Genomeditech) (left Y-axis, relative luminescence units), with inhibition percentages shown on the right Y-axis. Data are shown by drug mass concentration.

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
Fig 5. Response to Human APRIL / TNFSF13 Protein. The H\_TACI Reporter Cell Line (Cat. GM-C35041) at a concentration of 1E5 cells/well (96-well format) was stimulated with serial dilutions of Human APRIL/TNFSF13 Protein (ACRO/APL-H52D1) in assay buffer (RPMI 1640 + 1% FBS + 1% P.S) for 6 hours. The firefly luciferase activity was measured using the Luciferase Reporter Assay Kit (Genomeditech). The maximum induction fold was approximately [4.6]. Data are shown by drug mass concentration.

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
Fig 6. Response to Membrane Bound H\_APRIL(Trimer) HEK-293 Cell Line. H\_TACI Reporter Cell Line (Cat. GM-C35041) at a concentration of 1E5 cells/well (96-well format) was stimulated with serial dilutions of Membrane Bound H\_APRIL(Trimer) HEK-293 Cell Line (Cat. GM-C40503) for 6 hours. The firefly luciferase activity was measured using the Luciferase Reporter Assay Kit (Genomeditech).

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
Fig 7. The passage stability of response to Human BAFF Protein; His Tag. The passage 7, 17, and 27 of H\_TACI Reporter Cell Line (Cat. GM-C35041) at a concentration of 1E5 cells/well (96-well format) were stimulated with serial dilutions of Human BAFF Protein (Cat. GM-87735RP) in assay buffer (RPMI 1640+1% FBS+1% P.S) for 6 hours. The firefly luciferase activity was measured using the Luciferase Reporter Assay Kit (Genomeditech).

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
Fig 8. The passage stability of response to Anti-BAFF hlgG1 Antibody(belimumab). Serial dilutions of Anti-BAFF hlgG1 Antibody(belimumab) (Cat. GM-87690AB) was incubated with 12 ng/well of Human BAFF Protein; His Tag (Cat. GM-87735RP) for 1 hour in assay buffer (RPMI 1640+1% FBS+1% P.S). Subsequently, the passage 11, 21, and 31 of H\_TACI Reporter Cell Line (Cat. GM-C35041) at a concentration of 1E5 cells/well was added, and the coculture proceeded for an additional 6 hours. Firefly luciferase activity is then measured using the Luciferase Reporter Assay Kit (Genomeditech). Data are shown by drug mass concentration.

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
Fig 9. The passage stability of the H\_TACI Reporter Cell Line (Cat. GM-C35041) was determined by flow cytometry using APC anti-human CD267 (TACI) Antibody (Biolegend/311911).

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CD40: CD40L	
H_CD40(TNFRSF5) Reporter 293 Cell Line	H_CD40(TNFRSF5) Reporter Jurkat Cell Line
Cynomolgus_CD40 CHO-K1 Cell Line	Cynomolgus_CD40L CHO-K1 Cell Line
H_CD40(TNFRSF5) CHO-K1 Cell Line	H_CD40(TNFRSF5) HEK-293 Cell Line
H_CD40L CHO-K1 Cell Line	H_CD40L HEK-293 Cell Line
Anti-H_CD40 hIgG1 Antibody(APX005M)	Anti-H_CD40 hIgG1 Antibody(ravagalimab)
Anti-H_CD40L hIgG1 Antibody(dapirolizumab)	Anti-H_CD40L hIgG1 Antibody(frexalimab)
Biotinylated Human CD40 Protein; His-Avi Tag	Cynomolgus CD40 Protein; His Tag
Human CD40 Protein; His Tag	Human CD40L Protein; His Tag
IFN- $\alpha$	
IFN $\alpha$ Reporter HEK-293 Cell Line	IFN $\alpha$ Reporter MDCK Cell Line
IFN $\alpha$ Reporter THP1 Cell Line	
BCMA:BAFFR:TACI	
H_BCMA HEK-293 Cell Line	H_BAFFR Reporter Cell Line
H_BCMA Reporter Cell Line	Cynomolgus_BCMA CHO-K1 Cell Line
H_BCMA CHO-K1 Cell Line	
Anti-BAFF hIgG1 Antibody(belimumab)	Anti-BAFFR hIgG1 Antibody(ianalumab)
Anti-BCMA hIgG1 Antibody(Belantamab)	Anti-BCMA hIgG1 Antibody(SEA-BCMA)
Anti-BCMA hIgG4 Antibody(BCMB69)	
Biotinylated Human BAFF Protein; His-Avi Tag	Cynomolgus BAFF Protein; His Tag
Human BAFF Protein; His Tag	Mouse BAFF Protein; His Tag
BDCA2(CLEC4C)	
H_BDCA2 Reporter Jurkat Cell Line	Cynomolgus_BDCA2 CHO-K1 Cell Line
Cynomolgus_BDCA2 Jurkat Cell Line	H_BDCA2 CHO-K1 Cell Line
H_BDCA2 HEK-293 Cell Line	H_BDCA2 Jurkat Cell Line
Anti-H_BDCA2 hIgG1 Antibody(Litifilimab)	
Cynomolgus BDCA2 Protein; His Tag	Human BDCA2 Protein; His Tag
CD3	
Jurkat CD3-BsAb Reporter Cell Line	Cynomolgus_CD3 HEK-293 Cell Line
Cynomolgus_CD3E(Membrane Bound ECD) CHO-K1 Cell Line	H_CD3 CHO-K1 Cell Line
H_CD3 HEK-293 Cell Line	H_CD3E(Membrane Bound ECD) CHO-K1 Cell Line
Mouse_CD3 HEK-293 Cell Line	
Anti-CD3 epsilon hIgG1 Antibody [OKT-3 (muromonab)]	Anti-CD3 hIgG1 Antibody(CH2527)

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