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H_VEGF Reporter 293 Cell line

H_VEGF Reporter 293 细胞

For research use only!

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Reagent	Specification	Manufacturer/Catalogue No.
G418	10 mg	Genomeditech/GM-040402-1
Blasticidin	10 mg	Genomeditech/GM-040404-1
Pen/Strep	100 mL	Thermo/15140-122
Fetal Bovine Serum	500 mL	ExCell/FSP500
DMEM	500 mL	gibco/C11995500BT
96 Well Clear V-Bottom Tissue Culture	96-well	Corning/3894
96 well round well culture plate	96-well	NEST/701001
96 well White Flat Bottom Polystyrene Not Treated Microplate	96-well	Corning/3912
Human VEGF 121	10 µg	Novoprotein/C744
Human VEGF165	10 µg	Novoprotein/C083
Human KDR(VEGFR2)	10 µg	Novoprotein/CJ92
Anti-VEGF hIgG1 Antibody(Bevacizumab)	/	Genomeditech/GM-51978AB
Recombinant Mouse VEGF 164	10 µg	Novoprotein/CX73
Human VEGF165 Protein; His Tag	10 µg	Genomeditech/GM-88072RP
Anti-VEGFR2 hIgG1 Antibody(ramucirumab)	/	Genomeditech GM-80717AB
GMOne-Step 2.0 Luciferase Reporter Gene Assay Kit	/	Genomeditech/GM-040513

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Equipment	Manufacturer/Catalogue No.
细胞计数仪	ThermoFisher Scientific/Countess 3
酶标仪	Moleculardevices/SpectraMax L

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
Fig 2. Response to Recombinant Human VEGF 121. The H_VEGF Reporter 293 Cell Line (Cat. GM-C09057) at a concentration of 1.5E4 cells/well (96-well format) was stimulated with serial dilutions of Recombinant Human VEGF 121 (Novoprotein/C744) in assay buffer (DMEM + 1% FBS + 1% P.S) for 7 hours. The firefly luciferase activity was measured using the Luciferase Reporter Assay Kit (Genomeditech). The maximum induction fold was approximately [67.2]. Data are shown by drug mass concentration.

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
Fig 3. Response to Recombinant Human VEGF165. The H_VEGF Reporter 293 Cell Line (Cat. GM-C09057) at a concentration of 1.5E4 cells/well (96-well format) was stimulated with serial dilutions of Recombinant Human VEGF165 (Novoprotein/C083) in assay buffer (DMEM + 1% FBS + 1% P.S) for 7 hours. The firefly luciferase activity was measured using the Luciferase Reporter Assay Kit (Genomeditech). The maximum induction fold was approximately [9.5]. Data are shown by drug mass concentration.

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
Fig 4. Response to Anti-VEGF hIgG1 Antibody(Bevacizumab). Serial dilutions of Anti-VEGF hIgG1 Antibody(Bevacizumab) (Cat. GM-51978AB) was incubated with 0.5 ng/well of Recombinant Human VEGF165 (Novoprotein/C083) for 1 hour in assay buffer (DMEM + 1% FBS + 1% P.S). After pre-incubation, add the mixture to the H_VEGF Reporter 293 Cell Line (Cat. GM-C09057) at a density of 1.5E4 cells/well in a 96-well format, and incubate for 6 hours. The firefly luciferase activity was measured using the Luciferase Reporter Assay Kit (Genomeditech). The results indicated maximum blocking folds of approximately [6.8]. Data are shown by drug mass concentration.

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
Fig 5. Response to Recombinant Human VEGFR2 (C-Fc). Serial dilutions of Recombinant Human VEGFR2 (C-Fc) (Novoprotein/CJ92) was incubated with 1 ng/well of Recombinant Human VEGF165 (Novoprotein/C083) for 1 hour in assay buffer (DMEM + 1% FBS + 1% P.S). After pre-incubation, add the mixture to the H_VEGF Reporter 293 Cell Line (Cat. GM-C09057) at a density of 1.5E4 cells/well in a 96-well format, and incubate for 7 hours. The firefly luciferase activity was measured using the Luciferase Reporter Assay Kit (Genomeditech). The results indicated maximum blocking folds of approximately [16.2]. Data are shown by drug molar concentration.

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
Fig 6. The passage stability of response to Recombinant Human VEGF 121. The passage 5, 12, 22 and 32 of H_VEGF Reporter 293 Cell Line (Cat. GM-C09057) at a concentration of 1.5E4 cells/well (96-well format) were stimulated with serial dilutions of Recombinant Human VEGF 121 (Novoprotein/C744) in assay buffer (DMEM + 1% FBS + 1% P.S) for 7 hours. The firefly luciferase activity was measured using the Luciferase Reporter Assay Kit (Genomeditech). Data are shown by drug molar concentration.

附录 3: Mouse VEGF164 验证


Fig 7. Response to Recombinant Mouse VEGF164. The H_VEGF Reporter 293 Cell Line (Cat. GM-C09057) at a concentration of 1.5E4 cells/well (96-well format) was stimulated with serial dilutions of Recombinant Mouse VEGF164 (Novoprotein/CX73) in assay buffer (DMEM + 1% FBS + 1% P.S) for 6 hours. The firefly luciferase activity was measured using the Luciferase Reporter Assay Kit (Genomeditech). The maximum induction fold was approximately [16.5]. Data are shown by drug mass concentration.

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
Fig 8. H_VEGF Reporter 293 Cell Line (Cat. GM-C09057) was determined by flow cytometry using Anti-VEGFR2 hIgG1 Antibody(ramucirumab) (Cat. GM-80717AB).

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
Fig 9. Response to Human VEGF165. The H_VEGF Reporter 293 Cell Line (Cat. GM-C09057) at a concentration of 1.5E4 cells/well (96-well format) was stimulated with serial dilutions of Human VEGF165 Protein; His Tag (Cat. GM-88072RP) in assay buffer (DMEM + 1% FBS + 1% P.S) for 7 hours. The firefly luciferase activity was measured using the Luciferase Reporter Assay Kit (Genomeditech). The maximum induction fold was approximately [17.1]. Data are shown by drug mass concentration.

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VEGF	
H_VEGF Reporter 293 DDX35TM Cell Line	
Anti-mouse VEGFR-2 mIgG2a Antibody(DC101)	Anti-mouse VEGFR-2 RIgG1 Antibody(DC101)
Anti-VEGF hIgG1 Antibody(Bevacizumab)	Anti-VEGF hIgG1 Reference Antibody (Bevbio)
Anti-VEGFxPD1 hIgG1 Reference Antibody (Ivobio)	

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