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H\_IL20 IL22 IL24 Reporter 293 Cell Line

H\_IL20 IL22 IL24 Reporter 293 细胞系

For research use only!

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Reagent	Specification	Manufacturer/Catalogue No.
DMEM	500 mL	gibco/C11995500BT
Fetal Bovine Serum	500 mL	ExCell/FSP500
Pen/Strep	100 mL	Thermo/15140-122
Blasticidin	10 mg	Genomeditech/GM-040404-1
G418	1 g	Genomeditech/GM-040402-1
Hygromycin	1 g	Genomeditech/GM-040403-1
Puromycin	25 mg	Genomeditech/GM-040401-1
96 Well Clear V-Bottom Tissue Culture	96-well	Corning/3894
96 well round well culture plate	96-well	NEST/701001
96 well White Flat Bottom Polystyrene Not Treated Microplate	96-well	Corning/3912
Human IL22 Protein	/	Genomeditech/GM-88080RP
Recombinant Human IL-20 Protein	20 $\mu\text{g}$	SinoBiological/13060-HNAE
Recombinant Human IL-24 Protein (His Tag)	50 $\mu\text{g}$	SinoBiological/12265-H08H
IL-19 Protein, Human	10 $\mu\text{g}$	UA BIOSCIENCE/UA040152
Anti-IL-22RA hIgG1 Reference Antibody (Temptobio)		Genomeditech/GM-88094MAB
Anti-H_IL-22 hIgG1 Antibody(Fezakinumab)		Genomeditech/GM-46509AB
PE anti-human CD210b (IL-10RB) Antibody		Biologend/396803
GMOne-Step 2.0 Luciferase Reporter Gene Assay Kit	1000 tests	Genomeditech/GM-040513C

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Equipment	Manufacturer/Catalogue No.
细胞计数仪	ThermoFisher Scientific/Countess 3
酶标仪	Moleculardevices/SpectraMax L

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
Fig 2. Response to human IL20, IL22, IL24 and IL19 protein. The H\_IL20 IL22 IL24 Reporter 293 Cell Line (Cat. GM-C41955) at a concentration of  $1.5 \times 10^4$  cells/well (96-well format) was stimulated with serial dilutions of Recombinant Human IL-20 Protein (SinoBiological/13060-HNAE), Human IL-22 Protein; His Tag (Cat. GM-88080RP), Recombinant Human IL-24 Protein (His Tag) (SinoBiological/12265-H08H) and IL-19 Protein, Human (UA BIOSCIENCE/UA040152) in assay buffer (DMEM + 1% FBS + 1% P.S) for 16 hours. The firefly luciferase activity was measured using the Luciferase Reporter Assay Kit (Genomeditech). The maximum induction folds of IL20, IL22, IL24 were approximately [1156.4], [901.3] and [1205.04]. Data are shown by drug mass concentration.

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
Fig 3. Inhibition of Human IL-22 Protein-induced reporter activity by Anti-H\_IL-22 hIgG1 Antibody(Fezakinumab). Serial dilutions of Anti-H\_IL-22 hIgG1 Antibody(Fezakinumab) (Cat. GM-46509AB) were incubated with 0.5 ng/well of Human IL-22 Protein; His Tag (Cat. GM-88080RP) for 1 hour in assay buffer (DMEM + 1% FBS + 1% P.S). After pre-incubation, add the mixture to the H\_IL20 IL22 IL24 Reporter 293 Cell Line (Cat. GM-C41955) at a density of 1.5E4 cells/well in a 96-well format, and incubate for 16 hours. Firefly luciferase activity was measured using the Luciferase Reporter Assay Kit (Genomeditech) (left Y-axis, relative luminescence units, RLU), with inhibition percentages shown on the right Y-axis. Data are shown by drug mass concentration.


Fig 4. Inhibition of Human IL-22 protein-induced reporter activity by Temtokibart. Serial dilutions of the Anti-IL-22RA hIgG1 Reference Antibody (Temtobio) (Cat. GM-88094MAB) were incubated with 1.5E4 cells/well of the H\_IL20 IL22 IL24 Reporter 293 Cell Line (Cat. GM-C41955) in a 96-well plate for 1 hour in assay buffer (DMEM + 1% FBS + 1% P.S). Subsequently, the Human IL-22 Protein; His Tag (Cat. GM-88080RP) at a density

of 0.5 ng/well was added, and the coculture proceeded for an additional 16 hours. Firefly luciferase activity was then measured using the Luciferase Reporter Assay Kit (Genomeditech) (left Y-axis, relative luminescence units), with inhibition percentages shown on the right Y-axis. Data are shown by drug mass concentration.


Fig 5. Inhibition of Human IL-20 protein-induced reporter activity by Temtokibart. Serial dilutions of the Anti-IL-22RA hIgG1 Reference Antibody (Temtobio) (Cat. GM-88094MAB) were incubated with 1.5E4 cells/well of the H\_IL20 IL22 IL24 Reporter 293 Cell Line (Cat. GM-C41955) in a 96-well plate for 1 hour in assay buffer (DMEM + 1% FBS + 1% P.S). Subsequently, the Recombinant Human IL-20 Protein (Sino Biological/13060-HNAE) at a density of 0.5 ng/well was added, and the coculture proceeded for an additional 16 hours. Firefly luciferase activity was then measured using the Luciferase Reporter Assay Kit (Genomeditech) (left Y-axis, relative luminescence units), with inhibition percentages shown on the right Y-axis. Data are shown by drug mass concentration.


Fig 6. Inhibition of Human IL-24 protein-induced reporter activity by Temtokibart. Serial dilutions of the Anti-IL-22RA hIgG1 Reference Antibody (Temtobio) (Cat. GM-88094MAB) were incubated with 1.5E4 cells/well of the H\_IL20 IL22 IL24 Reporter 293 Cell Line (Cat. GM-C41955) in a 96-well plate for 1 hour in assay buffer (DMEM + 1% FBS + 1% P.S). Subsequently, the Recombinant Human IL-24 Protein (His Tag) (SinoBiological/12265-H08H) at a density of 25 pg/well was added, and the coculture proceeded for an additional 16 hours. Firefly luciferase activity was then measured using the Luciferase Reporter Assay Kit (Genomeditech)

(left Y-axis, relative luminescence units), with inhibition percentages shown on the right Y-axis. Data are shown by drug mass concentration.

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
Fig 7. H\_IL20 IL22 IL24 Reporter 293 Cell Line (Cat. GM-C41955) was determined by flow cytometry using Anti-IL-22RA hIgG1 Reference Antibody (Temtobio) (Cat. GM-88094MAB).


Fig 8. H\_IL20 IL22 IL24 Reporter 293 Cell Line (Cat. GM-C41955) was determined by flow cytometry using PE anti-human CD210b (IL-10RB) Antibody (Biolegend/396803).

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
Fig 9. The mRNA expression levels of H\_IL20RB in the H\_IL20 IL22 IL24 Reporter 293 Cell Line (Cat. GM-C41955) were determined by RT-qPCR.

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IL-22	
H_IL22 Reporter 293 Cell Line	Cynomolgus_IL22R HEK-293 Cell Line
H_IL22R CHO-K1 Cell Line	H_IL22R HEK-293 Cell Line
Mouse_IL22R HEK-293 Cell Line	
Anti-H_IL-22 hIgG1 Antibody(Fezakinumab)	Anti-H_IL-22R1 hIgG1 Antibody(ARGX-112)
Biotinylated Human IL-22 Protein; His-Avi Tag	Human IL-22BP Protein; hFc Tag
Mouse IL-22RA1 Protein; His Tag	

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