

Cynomolgus B7H4 Protein; His Tag


Product Information

Product Name	Cynomolgus B7H4 Protein; His Tag
Storage temp.	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles.
Catalog# / Size	GM-84641RP-100 / 100 μg GM-84641RP-1000 / 1 mg

Protein Information

Alternative Names	B7-H4,VTCN1,B7S1,B7h.5,B7H4B7X,RO1291
Source	Cynomolgus B7H4 Protein; His Tag (GM-84641RP) is expressed from human 293 cells (HEK-293). It contains AA Phe29-Ala258 (Accession # F7B770-1). This protein carries a His tag at the C-terminus.
Purity	> 90% as determined by SDS-PAGE
Endotoxin	< 1 EU/ μg , determined by LAL gel clotting assay
Predicted Mol Mass	27.2 KDa
Formulation	Supplied as a 0.2 μm filtered solution of PBS, pH7.4.
Description	B7H4 protein is a cell surface protein belonging to the B7/CD28 immune regulatory family. This protein plays a crucial role in regulating immune responses and tumor immune evasion. Overexpression of B7H4 is associated with the development and progression of various tumors, making it a potential target for cancer therapy. Additionally, B7H4 may also promote tumor immune evasion by influencing the function of other immune cells in the tumor microenvironment, such as inhibiting natural killer cells (NK cells) and regulatory T cells (Treg cells). The B7H4 protein is of great interest in the fields of immunotherapy and tumor immunotherapy.

SDS-PAGE




On SDS-PAGE under reducing (R) condition. The gel was stained overnight with Coomassie Blue. The purity of the protein is greater than 90%.

Bioactivity-ELISA

Cynomolgus B7H4 Protein; His Tag ELISA

0.2 μ g Cynomolgus B7H4 Protein; His Tag of per well



Cynomolgus B7H4 Protein; His Tag (Catalog # GM-84641RP) was immobilized at 2 μ g/ml (100 μ L/well). Increasing concentrations of Anti-H_B7H4(VTCN1) hIgG1 Antibody (FPA150) (Catalog # GM-24028AB) were added.