

Anti-TREM1 hIgG1 Antibody(Clone A)

Product information

GM-88441AB-10	10 µg
GM-88441AB-100	100 µg
GM-88441AB-1000	1 mg

Antibody Information

Species Reactivity	Human; Cynomolgus
Clone	Clone A
Source/Isotype	Human IgG1, Kappa
Application	Bioactivity-ELISA
Target	Detects TREM1
Gene	TREM1
Other Names	CD354, TREM-1
Gene ID	54210 (Human);102129811 (Cynomolgus)
Background	<p>Targeting TREM1 (Triggering Receptor Expressed on Myeloid cells 1) with specific antibodies holds significant therapeutic promise. TREM1 is a key amplifier of innate immune responses, and its overactivation drives pathological inflammation in sepsis, ischemic stroke, inflammatory bowel disease, and COVID-19.</p> <p>Anti-TREM1 antibodies can selectively block this amplification loop without completely abolishing host defense, offering a safer immunomodulatory strategy than broad anti-inflammatory agents. Emerging evidence also implicates TREM1 in tumor immune evasion and fibrosis, expanding potential applications to oncology and chronic organ diseases. Thus, developing anti-TREM1 biologics represents a valuable approach for “resetting” dysregulated inflammation across diverse conditions.</p>
Storage	Store at 2-8°C short term (1-2 weeks). Store at ≤ -20°C long term. Avoid repeated freeze-thaw.
Formulation	Supplied as a 0.2 µm filtered solution of PBS, pH7.2-7.4.
Endotoxin	< 1 EU/mg, determined by LAL gel clotting assay

Data Examples

Bioactivity-ELISA

Human TREM1 Protein; His Tag (Catalog # GM-88422RP) was immobilized at 2 µg/ml (100 µL/well). Increasing concentrations of Anti-TREM1 hIgG1 Antibody(Clone A) (Catalog # GM-88441AB) were added.




Fig. ELISA

Bioactivity-ELISA

Cynomolgus TREM1 Protein; His Tag (Catalog # GM-88424RP) was immobilized at 2 µg/ml (100 µL/well). Increasing concentrations of Anti-TREM1 hIgG1 Antibody(Clone A) (Catalog # GM-88441AB) were added.




Fig. ELISA